

New Journal

```
Well = Screen.Status.WellName + " " + "Site" + str(Screen.Status.SiteNum)
Power = str(Component.405_Laser_Power.Position) + "% " + "Laser Power"
```

```
IF Screen.Status.WaveNum=1 THEN
  Segmentation of miRFP for Targeting
```

```
  EDIT ANALYSIS for MiRFP from here [LINE 1] to [Line 19]
```

```
  1: New "LowPass" = Basic Filters("Camera GFP", 5, 5)
  2: New "Segmentation" = Count Nuclei([1: Basic Filters])
  3: Close([1: Basic Filters])
  4: Threshold Image([2: Count Nuclei], 1, 255, Inclusive)
  5: New "miRFP Binary" = Binarize([2: Count Nuclei]), high = current value, low = current value
  6: Close([2: Count Nuclei])
  IF size_dilate = 0 THEN
```

```
    7: Create Regions Around Objects([5: Binary Operations])
    8: Transfer Regions([5: Binary Operations], "Camera GFP", ALLREGIONS)
    9: Close([5: Binary Operations])
```

```
  ELSE
```

```
    10: Integrated Morphometry - Load State("miRFP Segmentation_Single Bin")
    11: Integrated Morphometry - Measure([5: Binary Operations])
    12: Set Image Zoom([5: Binary Operations], 100)
    13: New "Centroid_Binary" = As Displayed([Last Result], Entire image)
    14: Set Color Threshold([13: As Displayed], RGBCOLORMODEL and INCLUSIVEHUE and STATEINCLUSIVE, 254,
    15: New "MiRFP Centroid" = Binarize([13: As Displayed]), 1, 255
    16: New "Dilate" = Morphological Dilate([15: Binary Operations], Circle)
        size = size_dilate
    17: Create Regions Around Objects("Dilate")
    18: Transfer Regions("Dilate", "Camera GFP", ALLREGIONS)
    19: Close([16: Morphological Dilate])
    20: Close([15: Binary Operations])
    21: Close([5: Binary Operations])
    22: Close([13: As Displayed])
```

```
  END IF
```

Make the ROI's the same Color for overlay purposes later (Example Hi Intensity vs. Low Signal)

```
23: Select Image("Camera GFP")
IF Image.NumRegions>=1 THEN
  FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1
    Region.ColorBlue = 255
    Region.ColorGreen = 0
    Region.ColorRed = 0
  NEXT
ELSE
  No Regions Present
END IF
```

\*\*\*\*\*TARGETING FOR MOSAIC -NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE\*\*\*\*\*

```
24: Select Image("Camera GFP")
IF Image.NumRegions >=1 THEN
  Reg = 1
  Total_Regions = Image.NumRegions +Total_Regions
  Regions are Saved and loaded to Mosaic for Targeting
  Run Journal C:\MM\app\mmproc\journals\Mosaic journals\Save ROI Batch 1.jnl
  current_illumination = Device.Illumination.Setting
  25: Select Illumination("Camera GFP")
  26: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region,
    Coordinate system setting = Device.Magnification.Setting
```

```
Mask Exposure Duration [ms] = pulsetime
27: Delay(MILLISEC)
Time = pulsetime
28: Select Illumination("Camera GFP")
Setting Name = current_Illumination
29: Clear All Regions("Camera GFP")
ELSE
Reg = 0
Dont' Utilize the Mosaic on this Field of View
END IF
ELSE
END IF
```

Create Overlay Images for Final Display Purposes

```
IF Screen.Status.WaveName="Camera RFP_532x_Post" THEN
IF Reg=1 THEN
30: New "Color Combine" = Color Combine("Camera RFP_532x_Post", "Camera GFP", [None])
31: Text([30: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 12, "%Well% Pulse Time: %pulsetime% Power: %Power%")
32: Load Regions([30: Color Combine], "Batch 1 ROI Targeting")
33: Add to "Target Regions Overlay" = As Displayed([30: Color Combine], Entire image)
34: Close([30: Color Combine])
ELSE
35: New "Color Combine" = Color Combine("Camera RFP_532x_Post", "Camera GFP", [None])
36: Text([35: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 12, "%Well% Pulse Time: %pulsetime% Power: %Power%")
37: Add to "Target Regions Overlay" = As Displayed([35: Color Combine], Entire image)
38: Close([35: Color Combine])
END IF
ELSE
END IF
```