

## 3 Level Segmentation-

Low Expression MiRFP in Blue, Middle in Green and High in Red-

Pulse1 is Correlated to Blue (Low) - 800ms recommended

Pulse2 Length is Correlated to Green (Middle) - 200ms recommended

Pulse3 Length is Correlated to Red (High) - 50ms recommended

Note that GFP/Dendra2 image must be taken after miRFP image. So, images should go GFP --> NIR --> RFP (post)

Well = Screen.Status.WellName + " " + "Site" + str(Screen.Status.SiteNum)

Power = str(Component.405\_Laser\_Power.Position) + "% " + "Laser Power"

IF Screen.Status.WaveName="Camera GFP" THEN

Reg1 = 0

Reg2 = 0

Reg3 = 0

Segmentation of miRFP for Targeting: Threshold level after segmentation to creat binay imageLow Gray (0-1) Middle Green (0-1) High Red (0-1)

1: New "LowPass" = Basic Filters("Camera NIR", 5, 5)

2: New "Segmentation" = Count Nuclei(Src=[1: Basic Filters])

3: Close([1: Basic Filters])

4: Threshold Image([2: Count Nuclei], 1, 65535, Inclusive)

5: New "miRFP Binary" = Binarize([2: Count Nuclei]), high = current value, low = current value

6: Close([2: Count Nuclei])

Segment Bottom Quartile

7: Integrated Morphometry - Load State("0-Q1")

8: Integrated Morphometry - Measure("Camera NIR", mask image = [5: Binary Operations])

9: Integrated Morphometry - Log Data("Camera NIR", OBJECTS, CURRENTDATA, 1, 2)

10: Integrated Morphometry - Create Objects Mask()

11: Create Regions Around Objects("IMA Objects Mask")

IF Image.NumRegions>=1 THEN

Reg1 = 1

FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1

Region.ColorBlue = 255

Region.ColorGreen = 0

Region.ColorRed = 0

NEXT

12: Save Regions("IMA Objects Mask", "Batch 1 ROI Targeting")

ELSE

No Regions Present at the Low Threshold

END IF

13: Close("IMA Objects Mask")

Segment Second Quartile Nuclei

14: Integrated Morphometry - Load State("Q1-Q2")

15: Integrated Morphometry - Measure("Camera NIR", mask image = [5: Binary Operations])

16: Integrated Morphometry - Log Data("Camera NIR", OBJECTS, CURRENTDATA, 1, 2)

17: Integrated Morphometry - Create Objects Mask()

18: Create Regions Around Objects("IMA Objects Mask")

IF Image.NumRegions>=1 THEN

Reg2 = 1

FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1

Region.ColorBlue = 0

Region.ColorGreen = 255

Region.ColorRed = 0

NEXT

19: Save Regions("IMA Objects Mask", "Batch 2 ROI Targeting")

ELSE

No Regions Present at the Low Threshold

END IF

20: Close("IMA Objects Mask")

Segment Third Quartile Nuclei

21: Integrated Morphometry - Load State("Q2-Q3")

22: Integrated Morphometry - Measure("Camera NIR", mask image = [5: Binary Operations])

23: Integrated Morphometry - Log Data("Camera NIR", OBJECTS, CURRENTDATA, 1, 2)

24: Integrated Morphometry - Create Objects Mask()

25: Create Regions Around Objects("IMA Objects Mask")

IF Image.NumRegions>=1 THEN

    Reg3 = 1

    FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1

        Region.ColorBlue = 0

        Region.ColorGreen = 0

        Region.ColorRed = 255

    NEXT

    26: Save Regions("IMA Objects Mask", "Batch 3 ROI Targeting")

ELSE

    No Regions Present at the Low Threshold

END IF

27: Close("IMA Objects Mask")

28: Close([5: Binary Operations])

END Segmentation of Image

\*\*\*\*\*TARGETING FOR MOSAIC -Batch 1 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

29: Select Image("Camera NIR")

IF Reg1=1 THEN

    Regions are Saved and loaded to Mosaic for Targeting

    30: Load Regions("Camera NIR", "Batch 1 ROI Targeting")

    current\_illumination = Device.Illumination.Setting

    31: Select Illumination("Camera GFP")

    32: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I

        Coordinate system setting = Device.Magnification.Setting

        Mask Exposure Duration [ms] = pulsetime

    33: Delay(MILLISEC)

        Time = pulsetime

    34: Select Illumination("Camera GFP")

        Setting Name = current\_illumination

    35: Clear All Regions("Camera NIR")

ELSE

    Reg1 = 0

    Dont' Utilize the Mosaic on this Field of View

END IF

36: Select Image("Camera NIR")

\*\*\*\*\*TARGETING FOR MOSAIC -Batch 2 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

IF Reg2=1 THEN

    Regions are Saved and loaded to Mosaic for Targeting

    37: Load Regions("Camera NIR", "Batch 2 ROI Targeting")

    current\_illumination = Device.Illumination.Setting

    38: Select Illumination("Camera GFP")

    39: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I

```

        Coordinate system setting = Device.Magnification.Setting
        Mask Exposure Duration [ms] = pulsetime2
40: Delay(MILLISEC)
    Time = pulsetime2
41: Select Illumination("Camera GFP")
    Setting Name = current_Illumination
42: Clear All Regions("Camera NIR")
ELSE
    Dont' Utilize the Mosaic on this Field of View
END IF

```

\*\*\*\*\*TARGETING FOR MOSAIC -Batch 3 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

```

43: Select Image("Camera NIR")
IF Reg3=1 THEN
    Regions are Saved and loaded to Mosaic for Targeting
44: Load Regions("Camera NIR", "Batch 3 ROI Targeting")
    current_illumination = Device.Illumination.Setting
45: Select Illumination("Camera GFP")
46: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I
        Coordinate system setting = Device.Magnification.Setting
        Mask Exposure Duration [ms] = pulsetime3
47: Delay(MILLISEC)
    Time = pulsetime3
48: Select Illumination("Camera GFP")
    Setting Name = current_Illumination
49: Clear All Regions("Camera NIR")
ELSE
    Reg1 = 0
    Dont' Utilize the Mosaic on this Field of View
END IF
50: Select Image("Camera NIR")

```

ELSE

END IF

Create Overlay Images for Final Display Purposes: Conditions: (1,1,1) (1,1,0) (1,0,1) (1,0,0) (0,1,1) (0,1,0) (0,0,1) (0,0,0)

```

IF Screen.Status.WaveName="Camera RFP_532x_Post" THEN
    IF Reg1=1 THEN
        IF Reg2=1 THEN
            IF Reg3=1 THEN
                Lo, Mid and High Regions Found Reg1=1 and Reg 2=1 and Reg3=1
                Reg1=1 and Reg2=1 and Reg3=1
                #1 (1,1,1)
                51: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
                52: Text([51: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well% Pulse Time Hi: %pulsetime% P
                53: Load Regions([51: Color Combine], "Batch 1 ROI Targeting")
                54: Load Regions([51: Color Combine], "Batch 2 ROI Targeting")
                55: Load Regions([51: Color Combine], "Batch 3 ROI Targeting")
                56: Add to "Target Regions Overlay" = As Displayed([51: Color Combine], Entire image)
                57: Close([51: Color Combine])
            ELSE
                #2(1,1,0)
                58: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
                59: Text([58: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well% Pulse Time Hi: %pulsetime% P

```

```

        60: Load Regions([58: Color Combine], "Batch 1 ROI Targeting")
        61: Load Regions([58: Color Combine], "Batch 2 ROI Targeting")
        62: Add to "Target Regions Overlay" = As Displayed([58: Color Combine], Entire image)
        63: Close([58: Color Combine])
    END IF
ELSE
    Reg1=1 and Reg 2=0 Reg3=1
    IF Reg3=1 THEN
        #3 (1,0,1)
        64: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        65: Text([64: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
        66: Load Regions([64: Color Combine], "Batch 1 ROI Targeting")
        67: Load Regions([64: Color Combine], "Batch 3 ROI Targeting")
        68: Add to "Target Regions Overlay" = As Displayed([64: Color Combine], Entire image)
        69: Close([64: Color Combine])
    ELSE
        #4 (1,0,0)
        70: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        71: Text([70: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
        72: Load Regions([70: Color Combine], "Batch 1 ROI Targeting")
        73: Add to "Target Regions Overlay" = As Displayed([70: Color Combine], Entire image)
        74: Close([70: Color Combine])
    END IF
END IF
ELSE
    Reg1=0 and Reg2=1
    IF Reg2=1 THEN
        Reg1=0 and Reg2=1 and Reg3=1
        IF Reg3=1 THEN
            #5 (0,1,1)
            75: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            76: Load Regions([75: Color Combine], "Batch 2 ROI Targeting")
            77: Load Regions([75: Color Combine], "Batch 3 ROI Targeting")
            78: Text([75: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
            79: Add to "Target Regions Overlay" = As Displayed([75: Color Combine], Entire image)
            80: Close([75: Color Combine])
        ELSE
            #6 (0,1,0)
            81: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            82: Load Regions([81: Color Combine], "Batch 2 ROI Targeting")
            83: Text([81: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
            84: Add to "Target Regions Overlay" = As Displayed([81: Color Combine], Entire image)
            85: Close([81: Color Combine])
        END IF
    ELSE
        IF Reg3=1 THEN
            #7 (0,0,1)
            86: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            87: Load Regions([86: Color Combine], "Batch 3 ROI Targeting")
            88: Text([86: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
            89: Add to "Target Regions Overlay" = As Displayed([86: Color Combine], Entire image)
            90: Close([86: Color Combine])
        ELSE
            #8 (0,0,0)
            91: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            92: Text([91: Color Combine], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %pulsetime%    Pu
            93: Add to "Target Regions Overlay" = As Displayed([91: Color Combine], Entire image)
            94: Close([91: Color Combine])
        END IF
    END IF
END IF

```

```
    END IF  
ELSE  
END IF
```