

3 Level Segmentation-

Variables and logs that must be pre-opened at the beginning of the experiment:

blankdatalog (data logging should be paused at beginning of experiment - note that image data will be saved to the datalog; first set of measurements are cytoplasm; second set are nucleus)

hi_value
pulsetime
md_value
pulsetime2
lo_value
pulsetime3

Pulse1 is Correlated to Blue (High) - 800ms recommended (intensely nuclear) - NEED TO DEFINE, variable = hi_value

Pulse2 Length is Correlated to Green (Middle) - 200ms recommended (medium nuclear) - NEED TO DEFINE, variable = md_value

Pulse3 Length is Correlated to Red (Low) - 50ms recommended (low nuclear) - NEED TO DEFINE, variable = lo_value

Note that GFP/Dendra2 image must be taken after mRFP image. So, images should go NIR --> GFP --> RFP (post)

Well = Screen.Status.WellName + "_" + "Site" + str(Screen.Status.SiteNum)

Power = str(Component.405_Laser_Power.Position) + "% " + "Laser Power"

IF Screen.Status.WaveName="Camera GFP" THEN

Crop GFP, maintaining original image name

name = "Camera GFP"+"_"+"Crop"

1: Create Region()

Position:X/From:X = 0

Position:Y/From:Y = 1

Width/To:X = Image.Width

Height/To:Y = Image.Height-2

2: Image/Plane("Camera GFP")

Image.Name = name

3: Clear All Regions([Current At Start])

Crop NIR, maintaining original image name

name = "Camera NIR"+"_"+"Crop"

4: Create Region()

Position:X/From:X = 0

Position:Y/From:Y = 1

Width/To:X = Image.Width

Height/To:Y = Image.Height-2

5: Image/Plane("Camera NIR")

Image.Name = name

6: Clear All Regions([Current At Start])

Begin segmentation

7: New "LowPass" = Basic Filters("Camera GFP_crop", 5, 5)

8: New "Segmentation" = Count Nuclei([7: Basic Filters])

9: Close([7: Basic Filters])

10: Threshold Image([8: Count Nuclei], 1, 65535, Inclusive)

11: New "mRFP Binary" = Binarize([8: Count Nuclei], high = current value, low = current value)

12: New "Dilate" = Dilate Image([11: Binary Operations], WITHOUTCLOSING, 1, 3, NODILATECLOSE)

Create image subsets containing only cytoplasm and nucleus pixels; these will be used to calculate stuff.

13: Create Regions Around Objects([12: Dilate Image])

14: New "Cytoplasm" = "Dilate" - [11: Binary Operations] + 0

15: New "NIR_and_cytoplasm" = "Camera NIR_crop" AND [14: Arithmetic]

IMPORTANT STEP: threshold the image such that pixels of 0 are not measured by the region analysis

16: Threshold Image([15: Arithmetic], 1, 65535, Inclusive)

17: New "Nucleus" = Erode Image([11: Binary Operations], 1, 1, KEEPLAST)

18: New "NIR_and_nucleus" = [17: Erode Image] AND "Camera NIR_crop"

IMPORTANT STEP: threshold the image such that pixels of 0 are not measured by the region analysis

19: Threshold Image([18: Arithmetic], 1, 65535, Inclusive)

Close un-needed steps (cytoplasm, nucleus, miRFP binary). KEEP Dilate (for regions) and NIR_and_cytoplasm/nucleus

20: Close([11: Binary Operations])

21: Close([14: Arithmetic])

22: Close([17: Erode Image])

Initialize region storage images - use thresholded segmentation image to give a template to shrink regions to fit nuclei

23: Image/Plane([8: Count Nuclei])

24: Rename Image([Last Result], "Reg1_storage")

Reg1 = 0

25: Image/Plane([8: Count Nuclei])

26: Rename Image([Last Result], "Reg2_storage")

Reg2 = 0

27: Rename Image([8: Count Nuclei], "Reg3_storage")

Reg3 = 0

Transfer regions from Dilate to cytoplasm_NIR so measurements can be made

28: Transfer Regions([12: Dilate Image], [15: Arithmetic], ALLREGIONS and CLEARSOURCE)

29: Close([12: Dilate Image])

30: Close("Camera GFP_crop")

31: Close("Camera NIR_crop")

32: Select Image([15: Arithmetic])

Loop through regions in NIR_and_cytoplasm: (1) measure intensitieis, ignoring unthresholded pixels (2) transfer to nucleus

IF Image.NumRegions > 0 THEN

33: Resume Data Logging()

FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1

34: Region Measurements([15: Arithmetic], CurrentPlane, Active Region, Milliseconds)

cytoplasm_avg_intensity = RegionMeasurements.Measurements.AverageIntensity

35: Transfer Regions([15: Arithmetic], [18: Arithmetic], SELECTEDREGION)

36: Region Measurements([18: Arithmetic], CurrentPlane, Active Region, Milliseconds)

nucleus_avg_intensity = RegionMeasurements.Measurements.AverageIntensity

First, ensure that cell has a readable miRFP value

Check if ratio is greater than hi value. If it is: transfer region to Segmentation; shrink region to fit; then transfer region to nucleus

IF cytoplasm_avg_intensity > 11000 THEN

IF nucleus_avg_intensity/cytoplasm_avg_intensity > hi_value THEN

37: Transfer Regions([18: Arithmetic], "Reg1_storage", SELECTEDREGION)

38: Shrink Region to Fit("Reg1_storage")

ELSE

Check if ratio is greater than md value. If it is, copy region to Reg2_storage. If not, check next threshold.

IF nucleus_avg_intensity/cytoplasm_avg_intensity > md_value THEN

39: Transfer Regions([18: Arithmetic], "Reg2_storage", SELECTEDREGION)

40: Shrink Region to Fit("Reg2_storage")

ELSE

Check if ratio is greater than hi value. If it is, copy region to Reg3_storage. If not, do not transfer region

IF nucleus_avg_intensity/cytoplasm_avg_intensity > lo_value THEN

41: Transfer Regions([18: Arithmetic], "Reg3_storage", SELECTEDREGION)

42: Shrink Region to Fit("Reg3_storage")

ELSE

END IF

END IF

END IF

ELSE

IF nucleus_avg_intensity > 11000 THEN

IF nucleus_avg_intensity/cytoplasm_avg_intensity > hi_value THEN

43: Transfer Regions([18: Arithmetic], "Reg1_storage", SELECTEDREGION)

44: Shrink Region to Fit("Reg1_storage")

ELSE

Check if ratio is greater than md value. If it is, copy region to Reg2_storage. If not, check next threshold

IF nucleus_avg_intensity/cytoplasm_avg_intensity > md_value THEN

45: Transfer Regions([18: Arithmetic], "Reg2_storage", SELECTEDREGION)

```

    46: Shrink Region to Fit("Reg2_storage")
  ELSE
    Check if ratio is greater than hi value. If it is, copy region to Reg3_storage. If not, do not transfer re
    IF nucleus_avg_intensity/cytoplasm_avg_intensity > lo_value THEN
      47: Transfer Regions([18: Arithmetic], "Reg3_storage", SELECTEDREGION)
      48: Shrink Region to Fit("Reg3_storage")
    ELSE
      END IF
    END IF
  END IF
ELSE
  END IF
END IF
NEXT
  49: Pause Data Logging()
ELSE
  END IF
50: Close([15: Arithmetic])
51: Close([18: Arithmetic])

Save regions in Reg1_storage, then close it
52: Select Image("Reg1_storage")
IF Image.NumRegions>=1 THEN
  Reg1 = 1
  FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1
    Region.ColorBlue = 255
    Region.ColorGreen = 0
    Region.ColorRed = 0
  NEXT
  53: Save Regions("Reg1_storage", "Batch 1 ROI Targeting")
ELSE
  No Regions Present at the Low Threshold
END IF
54: Close("Reg1_storage")

Save regions in Reg2_storage, then close it
55: Select Image("Reg2_storage")
IF Image.NumRegions>=1 THEN
  Reg2 = 1
  FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1
    Region.ColorBlue = 0
    Region.ColorGreen = 0
    Region.ColorRed = 255
  NEXT
  56: Save Regions("Reg2_storage", "Batch 2 ROI Targeting")
ELSE
  No Regions Present at the Low Threshold
END IF
57: Close("Reg2_storage")

Save regions in Reg3_storage, then close it
58: Select Image("Reg3_storage")
IF Image.NumRegions>=1 THEN
  Reg3 = 1
  FOR Image.ActiveRegion = 1 TO Image.NumRegions STEP 1
    Region.ColorBlue = 0
    Region.ColorGreen = 255
    Region.ColorRed = 0
  NEXT
  59: Save Regions("Reg3_storage", "Batch 3 ROI Targeting")

```

```
ELSE
    No Regions Present at the Low Threshold
END IF
60: Close("Reg3_storage")
END Segmentation of Image
```

*****TARGETING FOR MOSAIC -Batch 1 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

```
61: Select Image("Camera NIR")
IF Reg1=1 THEN
    Regions are Saved and loaded to Mosaic for Targeting
    62: Load Regions("Camera NIR", "Batch 1 ROI Targeting")
    current_illumination = Device.Illumination.Setting
    63: Select Illumination("Camera GFP")
    64: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I
        Coordinate system setting = Device.Magnification.Setting
        Mask Exposure Duration [ms] = pulsetime
    65: Delay(MILLISEC)
        Time = pulsetime
    66: Select Illumination("Camera GFP")
        Setting Name = current_illumination
    67: Clear All Regions("Camera NIR")
ELSE
    Reg1 = 0
    Dont' Utilize the Mosaic on this Field of View
END IF
```

*****TARGETING FOR MOSAIC -Batch 2 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

```
68: Select Image("Camera NIR")
IF Reg2=1 THEN
    Regions are Saved and loaded to Mosaic for Targeting
    69: Load Regions("Camera NIR", "Batch 2 ROI Targeting")
    current_illumination = Device.Illumination.Setting
    70: Select Illumination("Camera GFP")
    71: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I
        Coordinate system setting = Device.Magnification.Setting
        Mask Exposure Duration [ms] = pulsetime2
    72: Delay(MILLISEC)
        Time = pulsetime2
    73: Select Illumination("Camera GFP")
        Setting Name = current_illumination
    74: Clear All Regions("Camera NIR")
ELSE
    Dont' Utilize the Mosaic on this Field of View
END IF
```

*****TARGETING FOR MOSAIC -Batch 3 ROI's.....NO NEED TO EDIT SCRIPT/CHANGE BELOW THIS LINE

```
75: Select Image("Camera NIR")
IF Reg3=1 THEN
```

```

Regions are Saved and loaded to Mosaic for Targeting
76: Load Regions("Camera NIR", "Batch 3 ROI Targeting")
current_illumination = Device.Illumination.Setting
77: Select Illumination("Camera GFP")
78: Targeted Illumination = Targeted Illumination(Illum setting=Camera GFP, Coord setting=20X APO, Active region, I
    Coordinate system setting = Device.Magnification.Setting
    Mask Exposure Duration [ms] = pulsetime3
79: Delay(MILLISEC)
    Time = pulsetime3
80: Select Illumination("Camera GFP")
    Setting Name = current_illumination
81: Clear All Regions("Camera NIR")
ELSE
    Reg1 = 0
    Dont' Utilize the Mosaic on this Field of Vlew
END IF
82: Select Image("Camera NIR")

```

```
ELSE
```

```
END IF
```

Create Overlay Images for Final Display Purposes: Conditions: (1,1,1) (1,1,0) (1,0,1) (1,0,0) (0,1,1) (0,1,0) (0,0,1) (0,0,0)
 IF Screen.Status.WaveName="Camera RFP_532x_Post" THEN

```

  IF Reg1=1 THEN
    IF Reg2=1 THEN
      IF Reg3=1 THEN
        Lo, Mid and High Regions Found Reg1=1 and Reg 2=1 and Reg3=1
        Reg1=1 and Reg2=1 and Reg3=1
        #1 (1,1,1)
        (Disabled)83: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        (Disabled)84: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well% Pulse Time Hi: %pu
        (Disabled)85: Load Regions([Nonexistent command], "Batch 1 ROI Targeting")
        (Disabled)86: Load Regions([Nonexistent command], "Batch 2 ROI Targeting")
        (Disabled)87: Load Regions([Nonexistent command], "Batch 3 ROI Targeting")
        (Disabled)88: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
        (Disabled)89: Close([Nonexistent command])
      ELSE
        #2(1,1,0)
        (Disabled)90: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        (Disabled)91: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well% Pulse Time Hi: %pu
        (Disabled)92: Load Regions([Nonexistent command], "Batch 1 ROI Targeting")
        (Disabled)93: Load Regions([Nonexistent command], "Batch 2 ROI Targeting")
        (Disabled)94: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
        (Disabled)95: Close([Nonexistent command])
      END IF
    ELSE
      Reg1=1 and Reg 2=0 Reg3=1
      IF Reg3=1 THEN
        #3 (1,0,1)
        (Disabled)96: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        (Disabled)97: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well% Pulse Time Hi: %pu
        (Disabled)98: Load Regions([Nonexistent command], "Batch 1 ROI Targeting")
        (Disabled)99: Load Regions([Nonexistent command], "Batch 3 ROI Targeting")
        (Disabled)100: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
        (Disabled)101: Close([Nonexistent command])
      ELSE
        #4 (1,0,0)

```

```

        (Disabled)102: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
        (Disabled)103: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %p
        (Disabled)104: Load Regions([Nonexistent command], "Batch 1 ROI Targeting")
        (Disabled)105: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
        (Disabled)106: Close([Nonexistent command])
    END IF
END IF
ELSE
    Reg1=0 and Reg2=1
    IF Reg2=1 THEN
        Reg1=0 and Reg2=1 and Reg3=1
        IF Reg3=1 THEN
            #5 (0,1,1)
            (Disabled)107: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            (Disabled)108: Load Regions([Nonexistent command], "Batch 2 ROI Targeting")
            (Disabled)109: Load Regions([Nonexistent command], "Batch 3 ROI Targeting")
            (Disabled)110: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %p
            (Disabled)111: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
            (Disabled)112: Close([Nonexistent command])
        ELSE
            #6 (0,1,0)
            (Disabled)113: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            (Disabled)114: Load Regions([Nonexistent command], "Batch 2 ROI Targeting")
            (Disabled)115: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %p
            (Disabled)116: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
            (Disabled)117: Close([Nonexistent command])
        END IF
    ELSE
        IF Reg3=1 THEN
            #7 (0,0,1)
            (Disabled)118: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            (Disabled)119: Load Regions([Nonexistent command], "Batch 3 ROI Targeting")
            (Disabled)120: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %p
            (Disabled)121: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
            (Disabled)122: Close([Nonexistent command])
        ELSE
            #8 (0,0,0)
            (Disabled)123: New "Color Combine" = Color Combine("Camera NIR", "Camera GFP", [None])
            (Disabled)124: Text([Nonexistent command], 10, 10, 0, 255, 0, Arial, Bold, 11, "%Well%    Pulse Time Hi: %p
            (Disabled)125: Add to "Target Regions Overlay" = As Displayed([Nonexistent command], Entire image)
            (Disabled)126: Close([Nonexistent command])
        END IF
    END IF
END IF
ELSE
END IF
END IF

```