

Startup Journal for NLS Experiments with Script v1.0.

Variables to assign: pulsetime; pulsetime2; pulsetime3; title of log journal; hi_value; md_value; lo_value

1: Delay(65, MINUTES)

2: pulsetime = 800

3: pulsetime2 = 200

4: pulsetime3 = 50

hi_value = 1.244

md_value = 1.079

lo_value = 0.964

Other variables (DO NOT CHANGE)

5: size_dilate = 0

6: Component.405_Laser_Power.Position = 100

Total_Regions = 0

7: Select Focus Device(Device Name = Leica Z Motor)

8: Start Continuous Focusing(0)